
	Culture Protocol for skin fibroblasts XP6BE cell line	 Version 01 – 04/27/2025 Author : J. Haidar Validation : I. Lagroye Page 1 / 4
--	--	---

Cell Culture

The SV40-immortalized skin fibroblast cell line XP6BE derived from a 19-year-old female patient with xeroderma pigmentosum (complementation group D), was obtained from the Coriell Institute (Camden, NJ, USA).




Cells are cultured in Dulbecco's Modified Eagle Medium/Nutrient Mixture F-12 with high glucose (DMEM-F12 ; Sigma-Aldrich, D6429) supplemented with 10% fetal bovine serum (FBS), 100 U/mL penicillin, and 100 µg/mL streptomycin, and maintained at 37°C in a humidified atmosphere containing 5% CO₂.

Passaging Protocol

1. **Prepare Materials**
 - Warm culture medium and Trypsin-EDTA to 37°C.
2. **Remove Old Medium**
3. **PBS Wash**
 - Rinse cells once with D-PBS (no Ca²⁺, Mg²⁺) (2 mL per for T75 Flask)
4. **Cell Detachment**
 - Add 1 mL of Trypsin-EDTA per T25 flask (adjust for other dish sizes).
 - Incubate for 2–5 minutes at 37°C, checking for detachment.
5. **Stop Reaction**
 - Add equal volume of pre-warmed medium to stop the trypsin reaction.
 - Gently pipette up and down to break clumps.
 - Count cells.
6. **Re-seeding**

Thawing Protocol

1. **Preparation**
 - Prewarm complete DMEM/F-12 medium (containing 10% fetal bovine serum [FBS] and antibiotics) and additional FBS to 37°C in a water bath.
2. **Thawing**
 - Retrieve frozen XP6BE cell vials from –80°C and immediately proceed to thawing.
 - In a sterile biosafety cabinet, quickly add ~1 mL of the prewarmed medium directly into the vial to initiate thawing while gently pipetting up and down to dilute DMSO and dislodge the cell pellet.

 	Culture Protocol for skin fibroblasts XP6BE cell line	 Version 01 – 04/27/2025 Author : J. Haidar Validation : I. Lagroye Page 2 / 4
--	--	---

- Transfer the cell suspension promptly into a 15 mL sterile conical tube to minimize exposure time to DMSO.

3. Dilution and Washing

- Slowly add 9 mL of prewarmed complete medium to the tube to dilute the residual DMSO.
- Centrifuge at $300 \times g$ (approximately 1200 rpm) for 5 minutes at room temperature.
- Carefully aspirate and discard the supernatant without disturbing the pellet.
- Resuspend the cell pellet in 4.5 mL of complete DMEM (with 10% FBS).
- Transfer the cell suspension to a sterile T25 culture flask.

4. FBS Supplementation

- Supplement the culture with 500 μ L of prewarmed FBS to achieve a final concentration of 20% FBS.
- Incubate the flask at 37°C in a humidified 5% CO₂ incubator.

Freezing Protocol

One confluent T75 flask yields approximately 4 cryovials.

1. Preparation




- Aspirate and discard the culture medium from a confluent T75 flask.
- Rinse cells gently with 5 mL of sterile PBS to remove residual serum. Discard the PBS.

2. Cell Detachment

- Add 1 mL of 0.05% trypsin-EDTA solution to the flask.
- Incubate for 5 minutes at 37°C to allow cell detachment.

3. Neutralization and Collection

- Add 9 mL of complete culture medium (DMEM/F-12 supplemented with 10% FBS and 1% penicillin-streptomycin) to neutralize the trypsin.
- Disperse cell clumps thoroughly by pipetting up and down 5–6 times.
- Transfer the entire cell suspension into a 15 mL sterile conical tube.

 	Culture Protocol for skin fibroblasts XP6BE cell line	 Version 01 – 04/27/2025 Author : J. Haidar Validation : I. Lagroye Page 3 / 4
--	--	---

4. Centrifugation and Washing

- Centrifuge at $300 \times g$ (approximately 1200 rpm) for 5 minutes.
- Discard the supernatant and gently resuspend the cell pellet in sterile FBS
- Repeat centrifugation at 1200 rpm for 5 minutes to wash.

5. Freezing Medium Preparation

- Discard the supernatant and resuspend the pellet in 4 mL of medium consisting of 90% FBS and 10% DMSO.
- Gently pipette to achieve a uniform cell suspension.

6. Aliquoting and Labeling

- Dispense 1 mL of the cell suspension into each sterile cryovial.
- Label each vial with cell type, passage number, and date of freezing.

7. Freezing and Storage

- Place the cryovials in a freezer container containing isopropanol and store at -80°C for 72 hours.
- After 3 days, transfer vials to long-term at -80°C freezer.

Transient transfections

Transient transfections are performed in T75 flasks using linear polyethyleneimine (PEI, MW 25,000; Polysciences, Inc., 23966) at a DNA:PEI ratio of 1:5.

A total of 15 μg of DNA is used per flask, consisting of 2.25 μg of the probe construct and 12.75 μg of pcDNA3.1(+) empty vector. Following overnight incubation, cells are detached, resuspended in phenol red-free DMEM/F-12 (Thermo Fisher Scientific, 21063-029), and seed into white 96-well plates with clear bottoms (Greiner Bio-One, Courtaboeuf, France) at a density of 40,000 cells/well (200 μL /well of a 2×10^5 cells/mL suspension).



**Culture Protocol for skin fibroblasts
XP6BE cell line**



Version 01 – 04/27/2025

Author : J. Haidar

Validation : I. Lagroye

Page 4 / 4