

Materials and Reagents

- **Rat Primary Cortical Astrocytes** (stored in liquid nitrogen)
- **Astrocyte Growth Medium**
 - 85% Dulbecco's Modified Eagle Medium (DMEM) with 4.5 g/L glucose
 - 15% Fetal Bovine Serum (FBS)
- **Dulbecco's Phosphate Buffered Saline (D-PBS)** (without calcium and magnesium)
- **Trypan Blue Stain** (for viability assessment)
- **Sterile 15 mL conical tubes**
- **Trypsin-EDTA (0.05%)**
- **Hemocytometer or Automated Cell Counter**
- **Tissue-culture treated flasks, plates, or Petri dishes (uncoated)**
- **37°C Water Bath**
- **37°C Incubator with 5% CO₂ and 90% Humidity**

Thawing Rat Primary Cortical Astrocytes

1. **Prepare Growth Medium:** Pre-warm Astrocyte Growth Medium to 37°C.
2. **Thaw Cells**
 - Remove the vial from liquid nitrogen.
 - Quickly thaw in a 37°C water bath, gently swirling the vial.
 - Do not submerge the cap; thaw within <2 minutes.
3. **Transfer Cells**
 - Transfer the cells into a 15 mL tube.
 - Rinse the vial with 1 mL of growth medium and add dropwise to the cells.
 - Slowly add 8 mL of growth medium, swirling gently.
4. **Centrifugation:** Spin at 250 × g for 5 min at room temperature.
5. **Resuspend Cells:** Aspirate the supernatant and resuspend in 2 mL of fresh growth medium.
6. **Cell Counting:** Determine cell viability using Trypan Blue Staining.
7. **Plating Cells:** Seed cells at 20,000 cells/cm² in a pre-warmed, uncoated culture dish.
8. **Incubation:** Place in a 37°C, 5% CO₂ incubator.
9. **Media Change:** Replace 50% of the medium every 4–5 days.
10. **Passage when 100% Confluent.**

Passaging Rat Primary Cortical Astrocytes

1. **Prepare Materials**
 - Warm Astrocyte Growth Medium and Trypsin-EDTA to 37°C.
2. **Remove Spent Medium**

 	Culture Protocol for Rat Primary Cortical Astrocytes	
		Version 01 - 03/13/2025 Author : E.Pinnet Page 2 / 2

- Collect spent medium in a **sterile tube** for later use.
- 3. **PBS Wash**
 - Rinse cells once with **D-PBS (no Ca²⁺, Mg²⁺)** (2 mL per **10 cm²** surface).
- 4. **Cell Detachment**
 - Add **3 mL of Trypsin-EDTA per T75 flask** (adjust for other dish sizes).
 - Incubate for **up to 5 minutes at 37°C**, checking for detachment.
- 5. **Stop Reaction**
 - Add **equal volume of collected medium** to stop the trypsin reaction.
 - Gently pipette up and down to break clumps.
- 6. **Centrifugation**
 - Transfer cells to a **15 mL tube**, spin at **250 × g for 5 min.**
- 7. **Resuspend Cells**
 - Use **pre-warmed growth medium**, mix gently, and count cells.
- 8. **Re-seeding**
 - Plate at **20,000 cells/cm²** in **uncoated tissue-culture flasks**.
 - Change **medium every 4–5 days**.

Growth Conditions

- **Temperature:** 37°C
 - **CO₂ Level:** 5%
 - **Humidity:** 90%
-